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
¹ Under certain classifications of risk, the World Health Organization (WHO), or local regulatory authorities recommend special Creutzfeldt-Jakob disease (CJD/TSE) inactivation processing procedures. Consult WHO and local regulations for further information.

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- 遵循以下標準：AAMI TIR12, AAMI TIR30, AAMI TIR34, ANSI/AAMI ST79, ANSI/AAMI ST81, ISO 17664 和 "Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling Guidance for Food and Drug Administration Staff".
- 清洗劑：3M Rapid Multi-Enzyme Cleaner(70503), Prolystica 2X Concentrate Enzymatic Presoak And Cleaner (1C33), Prolystica 2X Concentrate Neutral Detergent(1C32)
- 自動清洗器– Miele PG 8535

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1. AAMI TIR12 : Designing, testing, and labeling reusable medical devices for re-processing in health care facilities: A guide for medical device manufacturers.
2. AAMI TIR30 : A compendium of processes, materials, test methods, and acceptance criteria for cleaning reusable medical devices.
3. AAMI TIR34: Water for the reprocessing of medical devices
4. ANSI/AAMI ST77 : Containment devices for reusable medical device sterilization.
5. ANSI/AAMI ST79 : Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
6. ANSI/AAMI ST81 : Sterilization of medical devices—Information to be provided by the manufacturer for the processing of resterilizable medical devices.
7. ASTM F565 : Standard Practice for Care and Handling of Orthopedic Implants and Instruments.
8. ISO 11607-1: Packaging for terminally sterilized re-usable instruments .
9. ISO 15223-1 : Medical devices- Symbols to be used with medical device labels, labeling and information to be supplied - Part 1:General requirements.
10. ISO 17664: Sterilization of re-usable instruments -Information to be provided by the manufacturer for the processing of resterilizable re-usable instruments.
11. ISO 17665-1: Sterilization of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices.
12. Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling Guidance for Food and Drug Administration Staff, 2015.

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Please note that this Surgical Technique Guide has been authored in the English language. Any translations into other languages have not been reviewed or approved by United Orthopedic Corporation and their accuracy cannot be confirmed. Any translated guide should be reviewed carefully prior to use and questions regarding a Surgical Technique Guide should be directed to United Orthopedic Corporation at unitedorthopedic.com/contact

The CE mark is valid only if it is also printed on the product label.

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mdi Europa GmbH